



Sample name: **pyridoxine HCl**  
 Assay name: **pH-metric pKa**  
 Assay ID: **17I-11015**  
 Filename: **C:\Sirius\_T3\Mehtap\20170911\_test\_comp\_pKa\17I-11015\_pyridoxine HCl\_pH-metric pKa.t3r**

Experiment start time: **9/11/2017 4:00:50 PM**  
 Analyst: **Dorothy Levorse**  
 Instrument ID: **T311053**





## Sample

pyridoxine HCl concentration factor 0.964  
 pyridoxine HCl stoichiometry 1.000  
 Chloride stoichiometry 1.000  
 Base pKa 1 4.82  
 Acid pKa 2 8.90

## Sample graphs






## Other graphs



Sample name: **pyridoxine HCl**  
 Assay name: **pH-metric pKa**  
 Assay ID: **17I-11015**  
 Filename: **C:\Sirius\_T3\Mehtap\20170911\_test\_comp\_pKa\17I-11015\_pyridoxine HCl\_pH-metric pKa.t3r**

Experiment start time: **9/11/2017 4:00:50 PM**  
 Analyst: **Dorothy Levorse**  
 Instrument ID: **T311053**

## Other graphs (continued)



## Assay Model

### Settings

Settings	Value	Date/Time changed	Imported from
Sample name	pyridoxine HCl	9/11/2017 2:25:57 PM	User entered value
Sample by	Weight		Default value
Sample weight	0.003220 g	9/11/2017 2:51:35 PM	User entered value
Formula weight	205.64 g/mol	9/11/2017 2:25:57 PM	User entered value
Solubility	Unknown		Default value
Molecular weight	169.18	9/12/2017 11:33:17 AM	User entered value
Individual pKa ionic environments	No		Default value
Number of pKas	2	9/12/2017 11:40:37 AM	User entered value
Sample is a	Ampholyte	9/12/2017 11:35:04 AM	User entered value
pKa 1	2.69	9/12/2017 11:13:47 AM	C:\Sirius_T3\Mehtap-files\20170911_test_compounds_pKa
Type	Base		Default value
pKa 2	4.83	9/12/2017 11:35:02 AM	User entered value
Type	Acid	9/12/2017 11:35:14 AM	User entered value
logp (XH2 +)	-10.00		Default value
logP (neutral XH)	-10.00		Default value
logP (X -)	-10.00	9/11/2017 2:25:57 PM	User entered value
Stoichiometry	1.00000		Default value
Aprotic counterion name	Chloride		From standards.xml file
Stoichiometry	1.00	9/12/2017 11:30:18 AM	User entered value
Charge per counterion	-1		From standards.xml file

## Assay Settings

Setting	Value	Original Value	Date/Time changed	Imported from
<b>General Settings</b>				
Analyst name	Dorothy Levorse			
<b>Standard Experiment Settings</b>				
Number of titrations	1			



## Assay Settings

Sample name: **pyridoxine HCl**  
Assay name: **pH-metric pKa**  
Assay ID: **17I-11015**  
Filename: **C:\Sirius\_T3\Mehtap\20170911\_test\_comp\_pKa\17I-11015\_pyridoxine HCl\_pH-metric pKa.t3r**

Experiment start time: **9/11/2017 4:00:50 PM**  
Analyst: **Dorothy Levorse**  
Instrument ID: **T311053**

## Assay Settings (continued)

Setting	Value	Original Value	Date/Time changed	Imported from
Minimum pH	1.800			
Maximum pH	12.200			
pH step between points of	0.200			
Minimum titrant addition	0.00002 mL			
Maximum titrant addition	0.10000 mL			
Argon flow rate	100%			
<b>Advanced General Settings</b>				
Detect turbidity using	Spectrometer			
Monitor at a wavelength of	500.0 nm			
Absorbance threshold of	0.100			
Collect turbidity sensor data	No			
Stir after titrant addition for	5 seconds			
For titrant addition, stir at	15%			
Start titration using	Cautious pH adjust			
<b>Titrant Pre-Dose</b>				
Titrant pre-dose	None			
<b>Assay Medium</b>				
Cosolvent in use	No			
ISA water volume	1.50 mL			
Water added	Automatic			
After water addition, stir for	5 seconds			
At a speed of	10%			
<b>Sample Sonication</b>				
Sonicate	No			
<b>Sample Dissolution</b>				
Perform a dissolution stage	Yes			
Adjust and hold pH for dissolution	To start pH			
Stir to dissolve for	120 seconds			
For dissolution, stir at	10%			
<b>Carbonate purge</b>				
Perform a carbonate purge	No			
<b>Temperature Control</b>				
Wait for temperature	Yes			
Required start temperature	25.0°C			
Acceptable deviation	0.5°C			
Time to wait	60 seconds			
Stir speed of	10%			
<b>Titration 1</b>				
Titrate from	Low to high pH			
Adjust to start pH	Yes			
After pH adjust stir for	10 seconds			
<b>Data Point Stability</b>				
Stir during data point collection	Yes			
For point collection, stir at	15%			
Delay before data point collection	0 seconds			
Number of points to average	20 points			
Time interval between points	0.50 seconds			
Required maximum standard deviation	0.00500 dpH/dt			
Stability timeout after	60 seconds			
<b>Experiment cleanup</b>				
Adjust pH to cleanup	To start pH			
And then stir for	60 seconds			
For cleaning, stir at	25%			
Then add water volume	0.25 mL			
And then stir for	30 seconds			



## Assay Settings

Sample name: **pyridoxine HCl**  
Assay name: **pH-metric pKa**  
Assay ID: **17I-11015**  
Filename: **C:\Sirius\_T3\Mehtap\20170911\_test\_comp\_pKa\17I-11015\_pyridoxine HCl\_pH-metric pKa.t3r**

Experiment start time: **9/11/2017 4:00:50 PM**  
Analyst: **Dorothy Levorse**  
Instrument ID: **T311053**

## Calibration Settings

Setting	Value	Date/Time changed	Imported from
Four-Plus alpha	0.079	9/11/2017 4:00:50 PM	C:\Sirius_T3\HCl17111.t3r
Four-Plus S	1.0031	9/11/2017 4:00:50 PM	C:\Sirius_T3\HCl17111.t3r
Four-Plus jH	1.1	9/11/2017 4:00:50 PM	C:\Sirius_T3\HCl17111.t3r
Four-Plus jOH	-0.5	9/11/2017 4:00:50 PM	C:\Sirius_T3\HCl17111.t3r
Base concentration factor	1.015	9/12/2017 11:13:47 AM	C:\Sirius_T3\Mehtap-files\20170911_test_compounds_pKa\KOH17111.t3r
Acid concentration factor	1.002	9/11/2017 4:00:50 PM	C:\Sirius_T3\HCl17111.t3r

## Instrument Settings

Setting	Value	Batch Id	Install date
Instrument owner	Merck		
Instrument ID	T311053		
Instrument type	T3 Simulator		
Software version	1.1.3.0		
Dispenser module		T3DM1100253	3/31/2009 6:24:52 AM
Dispenser 0	Water		3/31/2009 6:25:05 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Water (0.15 M KCl)	8-18-17	9/8/2017 9:22:43 AM
Dispenser 2	Acid		3/31/2009 6:25:11 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Acid (0.5 M HCl)	166940	9/8/2017 9:21:27 AM
Dispenser 1	Base		3/31/2009 6:25:21 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Base (0.5 M KOH)	01/06/17	9/8/2017 9:20:03 AM
Dispenser 5	Cosolvent		3/31/2009 6:26:24 AM
Syringe volume	2.5 mL		
Firmware version	1.2.1(r2)		
Distribution valve 5	Distribution Valve		3/31/2009 6:28:19 AM
Firmware version	1.1.3		
Port A	Methanol (80%, 0.15 M KCl)	8-15-17	9/7/2017 3:40:52 PM
Dispenser 3	Buffer		8/3/2010 6:05:16 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Phosphate Buffer		8/15/2017 11:15:20 AM
Dispenser 6	Octanol		10/22/2010 11:52:43 AM
Syringe volume	0.5 mL		
Firmware version	1.2.1(r2)		
Titrant	Octanol		8/15/2017 11:16:41 AM
Titrator		T3TM1100153	3/31/2009 6:24:17 AM
Horizontal axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Vertical axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Chassis I/O firmware version	1.11 AI1DI0DO4 Norgren I/O		
Probe I/O firmware version	1.1.1		
Electrode	T3 Electrode	T3E0769	8/15/2017 10:21:54 AM
E0 calibration	-4.08 mV		9/11/2017 4:01:14 PM
Filling solution	3M KCl	KCL095	9/10/2017 10:38:19 AM
Liquids			
Wash 1	50% IPA:50% Water		9/11/2017 9:22:00 AM
Wash 2	0.5% Triton X-100 in H2O		9/11/2017 9:22:03 AM
Buffer position 1	pH7 Wash		9/11/2017 9:22:06 AM
Buffer position 2	pH 7		9/11/2017 9:22:10 AM
Storage position			9/11/2017 9:23:34 AM
Wash water	2.4e+003 mL	8-16-17	8/16/2017 4:10:31 PM
Waste	9.2e+003 mL		8/15/2017 10:25:02 AM
Temperature controller			8/5/2010 7:35:13 AM



## Assay Settings

Sample name: **pyridoxine HCl** Experiment start time: **9/11/2017 4:00:50 PM**  
Assay name: **pH-metric pKa** Analyst: **Dorothy Levorse**  
Assay ID: **17I-11015** Instrument ID: **T311053**  
Filename: **C:\Sirius\_T3\Mehtap\20170911\_test\_comp\_pKa\17I-11015\_pyridoxine HCl\_pH-metric pKa.t3r**

## Instrument Settings (continued)

Setting	Value	Batch Id	Install date
Turbidity detector			3/31/2009 6:24:45 AM
Spectrometer		072390	11/23/2010 12:22:28 PM
Dip probe		11086	
Wavelength coefficient A0	185.563		
Wavelength coefficient A1	2.17439		
Wavelength coefficient A2	-0.000285622		
Total lamp lit time	65:29:35		11/23/2010 12:22:28 PM
Calibrated on	9/6/2017 9:33:02 AM		
Integration time	11		
Scans averaged	10		
Autoloader		T3AL1100237	11/10/2015 10:34:13 AM
Left-right axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Front-back axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Vertical axis firmware version	1.17 AI1DI2DO2 Stepper 2		
Chassis I/O firmware version	1.11 AI1DI0DO4 Norgren I/O		
Configuration			
Alternate titration position	Titration position		
Alternate reference position	Reference position		
Maximum standard vial volume	3.50 mL		
Maximum alternate vial volume	25.00 mL		
Automatic action idle period	5 minute(s)		
Titrant tube volume	1.3 mL		
Syringe flush count	3.50		
Flowing wash pump volume	20.0 mL		
Flowing wash stir duration	5 s		
Flowing wash stir speed	30%		
Solvent wash stir duration	5 s		
Solvent wash stir speed	30%		
Surfactant wash stir duration	5 s		
Surfactant wash stir speed	30%		
E0 calibration minimum number of points	10		
E0 calibration maximum standard deviation	0.01500		
E0 calibration timeout period	60 s		
E0 calibration stir duration	5 s		
E0 calibration preparation stir speed	30%		
E0 calibration buffer wash stir duration	5 s		
E0 calibration buffer wash stir speed	30%		
E0 calibration reading stir speed	0%		
Spectrometer calibration stir duration	5 s		
Spectrometer calibration stir speed	30%		
Spectrometer calibration wash pump volume	20.0 mL		
Spectrometer calibration wash stir duration	5 s		
Spectrometer calibration wash stir speed	30%		
Overhead dispense height	10000		

## Refinement Settings

Setting	Value	Default value
Turbidity detection method	Spectrometer	Spectrometer
Turbidity wavelength to assess	500.0 nm	500.0 nm
Turbidity maximum absorbance	0.100	0.100
Turbidity probe threshold	50.00	50.00